



## REDGUM SELF PROPELLED QUARTZ WHEELCHAIR

**\$495.00 - \$510.00** Price range: \$495.00 through \$510.00



**SKU:** RG310A18-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Introducing the Redgum Self-Propelled Quartz Wheelchair – sleek design paired with optimal functionality. Weighing just 13.5kg, it offers a robust framework without compromising portability. The chair boasts comfortable padded armrests that lift, durable nylon upholstery available in two striking colours, and a removable seat cover for added convenience. With 24" solid tyres, removable wheels, and user-friendly swing-away leg rests this wheelchair promotes both mobility and comfort.

### FEATURES

- Durable yet lightweight aluminium frame in classic black
- Foldable frame and backrest for compact storage
- Elevating padded armrests for additional comfort
- Resilient nylon upholstery in either vibrant Blue or Red
- Detachable comfort seat cover for ease of cleaning
- Self-propulsion enabled with quick-release rear wheels; hand brakes included
- Adjustable and removable swing-away leg rests; featuring colour-coordinated calf straps
- Height-adjustable swing-up footrests to suit individual needs
- Convenient push handles for attendant propulsion
- Handy storage pocket located on the backrest

### SPECIFICATIONS

- **Dimensions:**
  - Overall: 700mm W x 1050mm D x 900mm H
  - Seat: 485mm W (19 inches) x 440mm D x 500mm H (from floor)

- Armrest: 240mm above seat
- Backrest: 380mm H
- Footrest: 370 - 470mm above seat
- **Wheels:** Ø200mm (front), Ø600mm (rear)
- **Product Weight:** 13.5kg
- **Max User Weight:** 115kg