



# REDGUM SIDE FOLDING TABLE - SUITS KATIE DAY CHAIR

\$85.00 Original price was: \$85.00.<u>\$77.00</u>Current price is: \$77.00.



SKU: RG63SFT

## PRODUCT DESCRIPTION

### **PURPOSE & BENEFITS**

A attachable and foldable side table designed to complement The Katie Day Chair High Back & Adjustable Chair.

## **FEATURES**

- Side folding table can be fitting to the:
  - Katie High Back Chair
  - Katie Low Back Chair
- Can be fitted on left or right side of Katie Chair.
- Ideal for holding small items(Drink bottle / Phone / Book / TV remote etc).
- Hard wearing laminated platform with plastic edging 400 x 200mm table platform.
- Table can be folded down when not in use.
- Easy to engage / disengage folding table.
- Large user friendly release lever.
- Strong lightweight construction.

### **SPECIFICATIONS**

Table Width: 200mmTable Depth: 400mm

• Maximum Weight Loading: 10kg (evenly distributed)

• Product Weight: 3kg