

## REDGUM SPACE SAVER CORNER SHOWER STOOL

~~\$85.00~~ Original price was: \$85.00. \$69.00 Current price is: \$69.00.



**SKU:** RG553



## PRODUCT DESCRIPTION

### Purpose & Benefits

The **Space Saver Corner Shower Stool** is designed to provide a secure and comfortable seat in compact shower spaces—particularly ideal for the elderly, expectant mothers, or anyone who finds standing for long periods challenging. Its easy no-tools assembly, angled reinforced seat, and adjustable height make it both safe and versatile. Lightweight and portable, yet robust, it offers stability and independence, while the corrosion-resistant build and lifetime warranty ensure lasting reliability.

### Features

- Easy to assemble – no tools required
- Patented clip-in leg attachment system for quick and secure setup
- Height adjustable via Sure Lock Adjustment—25 mm increments for personalised fit
- Improved leg tips enhance stability on wet surfaces
- Angled, reinforced plastic seat for extra comfort and durability
- Strong yet lightweight construction for ease of handling

- Easy-clean design suitable for wet environments
- Corrosion-resistant and supported by a lifetime warranty

## Specifications

- **Seat Material:** Polypropylene (plastic)
- **Leg Material:** Anodised aluminium
- **Leg Tips:** Flexible PVC
- **Overall Dimensions (Lowest Setting):** 430 mm (W) × 360 mm (D) × 490 mm (H)
- **Seat Width / Depth:** 395 mm / 295 mm
- **Adjustable Seat Height:** 490–615 mm (floor to top of seat), 25 mm increments
- **Product Weight:** 1.5 kg
- **Maximum User Weight:** 130 kg
- **Flat Packed Size:** 420 mm × 310 mm × 80 mm (L × W × H)