



REDGUM SPACE SAVER SHOWER CHAIR

\$119.00



SKU: RG555H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Designed in Australia the Space Saver Range is a durable and easily assembled solution for Shower and bathrooms.

The Space Saver Shower Chair offers quick tool free assembly and is space efficient for applications where a light weight and durable shower chair is required.

FEATURES

- Easy to assemble – No tools required.
Patented clip in leg attachment system.
- Adjustable in height with Sure Lock Adjustment.
- Rubber leg ends for stability and grip.
- Angled seat for maximum comfort.
- Reinforced plastic seat.
- Ribbed arm handles.
- Strong lightweight construction.
- Easy Clean Design.
- Lifetime warranty against corrosion.

SPECIFICATIONS

- Overall Dimensions (at lowest setting): 515mm(W) x 450mm (D) x825mm (H)

- Seat Depth: 355mm
- Seat Height: 485-610mm (from the floor to the top of the seat).
- Adjustable height in 25mm adjustment increments.
- Seat Handle: 155mm (From top of seat to top of handle).
- Product Weight: 3.1kg
- Maximum User Weight: 130KG
- Flat Packed Size (Box): 395mm(L) x 350mm (W) x 125mm (H)