




REDGUM SPACE SAVER SHOWER STOOL



SKU: RG554H2

PRODUCT DESCRIPTION

Purpose & Benefits of Redgum Space Saver Shower Stool

The Redgum Space Saver Shower Stool is designed to enhance safety and comfort in the bathroom for individuals with limited mobility, such as the elderly, those recovering from surgery, or individuals with disabilities. This compact stool fits well in smaller shower enclosures, providing a secure place to sit while bathing without occupying too much space.

Features

- **Compact Design:** Ideal for smaller shower spaces, this stool ensures that users have a comfortable and stable place to sit without crowding the area.
- **Adjustable Height:** Comes with adjustable legs that allow the stool to be set at the optimal height for the user, ensuring both comfort and safety.
- **Non-slip Rubber Feet:** Equipped with rubber feet that prevent the stool from slipping, even on wet surfaces, enhancing stability and safety.
- **Durable Construction:** Made from high-quality materials that are water-resistant and designed to withstand the humid conditions of a bathroom.
- **Easy to Clean:** The stool's non-porous surfaces can be easily wiped down, which helps maintain hygiene and cleanliness.
- **Compatible** with [RG53GP Gel Pad](#)

Specifications

- **Material:** PVC, Aluminium, Polypropylene
- **Dimensions:**
 - Overall Width: 515mm (at lowest setting)
 - Overall Depth: 450mm (at lowest setting)
 - Overall Height: 640mm (at lowest setting)
- **Safe Working Load:** 130kg

The Redgum Space Saver Shower Stool is an essential aid for anyone needing support in the shower, offering stability and comfort without taking up excessive space. Its adjustable and secure design makes it an excellent choice for enhancing bathroom safety and independence.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au