



REDGUM SPACE SAVER TOILET RAISER GEL PAD



SKU: RG51GP

PRODUCT DESCRIPTION

Purpose & Benefits of Redgum Space Saver Toilet Raiser Gel Pad

The Redgum Space Saver Toilet Raiser Gel Pad is designed to enhance comfort and support for individuals using toilet seat raisers. This gel pad is particularly useful for those with reduced mobility or who spend extended periods on the toilet, as it helps alleviate pressure points, reducing discomfort and the risk of pressure sores.

Features

- **Pressure Relief:** The gel pad is crafted with 10 mm thick gel to effectively distribute weight and relieve pressure, making prolonged sitting more comfortable.
- **Easy to Use:** Equipped with a washable adhesive base that ensures the pad stays securely in place on the toilet seat raiser.
- **Hygienic and Easy to Clean:** The gel pad's surface is easy to clean, supporting optimal hygiene in bathroom environments.
- **Compatibility:** Specifically designed to fit the [RG515 Space Saver Toilet Seat Raiser](#), ensuring a perfect fit and enhanced stability.

Specifications

- **Material:** Durable, hygienic gel material suitable for frequent use and cleaning.
- **Dimensions:** 370 mm (W) x 380 mm (D) x 10 mm (H), designed to cover the entire seating area of the toilet raiser.
- **Weight:** 1 kg

This gel pad is an excellent addition for anyone seeking additional comfort and support with their

existing toilet setup, especially suitable for the elderly or those with mobility impairments. It is designed to fit seamlessly with specific models of toilet seat raisers, ensuring enhanced user satisfaction and improved daily living conditions.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au