



REDGUM SPACE SAVER TOILET SEAT RAISER

\$129.00



SKU: RG515

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

A portable or permanent use Toilet Seat Raiser from the Redgum Space Saver range.

FEATURES

- Easy to assemble – No tools required.
- Patented clip in leg attachment system.
- Adjustable in height with Sure Lock Adjustment.
- Contoured toilet seat + cover flap.
- Removable splash guard.
- Leg tips for excellent stability.
- Ribbed arm handles.
- Strong lightweight construction.
- Easy clean design.
- Can be used with the [Space Saver Gel Pad](#)

SPECIFICATIONS

- Overall Dimensions (at lowest setting): 540mm (W) x 440mm (D) x 655mm (H)
- Seat Width: 470mm (between arm handles)
- Seat Depth: 400mm
- Seat Height: 500 – 625mm (from floor to top of seat), 25mm adjustment increments
- Seat Handle: 150mm (from top of seat to top of handle)
- Product Weight: 3.5kg
- Maximum User Weight: 130kg

