

REDGUM STANDARD BED RAIL

\$41.00



SKU: RG611



PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Assistive grab rail to aid user to get in and out of bed comfortably, securely and helping to maintain their independence.

Underbed Universal Bed Rail offers support when you have difficulty getting in and out of bed or adjusting your position once in bed.

Redgum Brand Removable Bed Rail is able to be used on all mattresses including lift beds.

FEATURES

- Easy to assemble - no tools required
- Easy to locate and remove for bed making
- Strong lightweight construction
- Epoxy powder coated mild steel frame
- Can pull apart for storage and transport

SPECIFICATIONS

- Length: 600mm
- Length: 550mm
- Length: 495mm
- Product Weight: Approx. 2kg