


## REDGUM TELESCOPIC RAMP PAIR

**\$285.00**



**SKU:** RAM76T



## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Introducing the Redgum Telescopic Ramp, designed with your convenience and safety in mind! This innovative ramp offers an easy-to-extend function, allowing it to reach up to 4 feet in length, providing effortless access for wheelchairs and mobility scooters.

Weighing just 4.0 kg, this ramp is lightweight yet robust, with a generous distributed weight capacity of 270 kg. Whether it's for personal use at home or professional applications, it can handle the demands of daily usage with ease.

One of the standout features of this ramp is its grooved non-slip surface. This specially designed texture ensures a firm grip for wheels, reducing the chance of slippage and promoting a safe and stable passage.

### FEATURES

- **Easy to Extend Function:** Allows for quick adjustment up to 4 feet.
- **Sold as a Pair:** Ensuring consistency and reliability.
- **Grooved Non-Slip Surface:** Enhances safety by providing a secure grip.
- **Lightweight Design:** At only 4.0 kg, it's easy to transport and store.
- **Robust Weight Capacity:** Supports a distributed weight of up to 270 kg.

### SPECIFICATIONS

	Dimension Name	Dimension
Length		1800mm

	Dimension Name	Dimension
Width		158mm
Weight		4.0 kg