


REDGUM TELESCOPIC RAMPS ECONOMY PAIR



SKU: RAM70

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing our Redgum Telescopic Ramps Economy Pair, engineered to assist with your mobility needs. These ramps are crafted with an Aluminum Frame and come equipped with a textured non-slip grip surface, ensuring safety and stability.

Sold as a pair, these ramps are built with quality and durability in mind. Each ramp is capable of supporting a distributed weight limit of up to 270 kg, making them suitable for a wide range of applications.

With an adjustable length that extends from a minimum of 1060mm to a maximum of 1800mm, and a width of 190mm, these ramps are designed to accommodate various needs. Each ramp weighs just 3.5kg, making them easy to transport and set up.

FEATURES

- **270KG Distributed Weight Limit:** Robust enough for heavy-duty applications.
- **Sold as a Pair:** Ensuring consistency and reliability.
- **Aluminium Frame:** Offers sturdiness and longevity.
- **Textured Non-Slip Grip Surface:** Adds an extra layer of safety.

SPECIFICATIONS

Part No	Min Length	Max Length	Width	Weight
Ram 70	1060mm	1800mm	190mm	3.5kg ea

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au