



REDGUM TOILET SURROUND FRAME ADJUSTABLE WIDTH

\$84.00 Original price was: \$84.00.<u>\$75.60</u>Current price is: \$75.60.



SKU: RG52TS

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Redgum Toilet Surround Frame Adjustable Width** provides safe and stable support for users when sitting or standing during toileting. Designed with dual grab handles and a lightweight, foldable frame, it is an ideal aid for individuals who require extra balance or mobility assistance in the bathroom.

FEATURES

- Excellent sitting and standing aid when toileting
- Dual grab handles assist with safe positioning and transfers
- Height adjustable legs with five settings for a custom fit
- Non-slip leg tips for added stability and safety
- Quick and easy tool-free setup
- Simple to position and remove for cleaning or use by others
- Foldable design for convenient storage or transport

SPECIFICATIONS

Specification	Detail
Width Overall (Unfolded)	560 mm
Width Between Handles	495 mm
Depth (Unfolded)	500 mm
Seat Heights (Floor to Top of Handle	e) 600 mm / 625 mm / 650 mm / 675 mm / 700 mm
Overall Folded	560 mm (L) x 90 mm (W) x 600 mm (H)



Specification Detail

Product Weight 3 kg
Maximum Supportive Load 130 kg