




REDGUM FOLDING TRANSIT WHEELCHAIR

\$440.00 - \$495.00 Price range: \$440.00 through \$495.00



SKU: N/A

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

This **Folding Deluxe Transit Wheelchair** has been specifically designed to be lightweight aluminium construction powder coated folding frame, making it easy to both store and transport. It has a number of features which accommodate this purpose. The transport chair has a compact design, with a half folding back and cranked push handles which assist to reduce storage space. It is available in a Flame crackle pattern red and blue in colour. It is also extremely easy to lift as a result of its lightweight nature.

This Folding Deluxe Transit Wheelchair has a number of features which make the chair more comfortable, having desk style padded armrests which can swing back or be removed for transferring of the occupant to either chair, bed etc. The Transport chair also has a padded backrest and seat belt for safety. This allows the chairs user

Having puncture proof solid fixed 30cm (12") rear wheels and also puncture proof front 8 inch front castors. The chair is suitable for both indoors and outdoor use and also provides its user with a high level of manoeuvrability. The chair's wheels, as well as its strong nylon canvas seat mean that the chair is durable, easy to care for and easy to clean.

FEATURES

- Lightweight Aluminium Construction Powder Coated Folding Frame
- Fold Down Back Rests
- Lift up Armrests
- Padded Arm Rests
- Removable Armrests

- 12"Rear PU Tyres
- 8"Front Rubber Tyres
- Hard Wearing Nylon Upholstery
- Swing Away Leg Rests
- Attendant Rear Control Hand brakes
- Manual Brake on the Rear Wheels
- Seat Belt Fitted
- Rear Storage Pocket on Backrest

SPECIFICATIONS

- Seat Width 43.18cm (17"); 48.26cm (19"); 56cm (22")
- Seat Depth 390mm
- Seat Height 460mm
- Max User Capacity 17/19" - 130kg
- Max User Capacity 22" - 180kg
- Wheelchair Weight 17/19" - 10.5kg
- Wheelchair Weight 22" - 13.0kg