



REDGUM VITTORIA HD DUAL MOTOR VIBRATION MASSAGE BLACK LEATHER



SKU: LC0205

PRODUCT DESCRIPTION

Purpose & Benefits

The **Redgum Vittoria HD Dual Motor Vibration Massage Lift Chair - Black Leather** provides exceptional comfort and support for users who require assistance standing or reclining. With dual motors for independent movement and built-in vibration massage, it combines style, functionality, and luxury for everyday relaxation and mobility support.

Features

- Lift, rest, and recline functionality for complete user comfort
- Dual motor system allows independent control of backrest and legrest
- Inbuilt vibration massage function for relaxation and improved circulation
- Luxurious real leather finish with PU back in black
- Pocket spring seat base for enhanced comfort and support
- Backlit handset for ease of use
- Fitted with high-quality German actuators for reliability
- Supports users up to 220 kg

Specifications

- **Overall Width:** 108 cm
- **Overall Depth:** 98 cm
- **Overall Height:** 108.5 cm
- **Seat Width:** 64.5 cm
- **Seat Depth:** 54 cm
- **Seat Height:** 50 cm
- **Product Weight:** 58 kg
- **Maximum Lifting Capacity:** 220 kg
- **Warranty:** Frame – 3 years, Electrics – 2 years, Leather – 2 years

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au