



### **REDGUM WALKING FRAME WITH FRONT WHEELS**

~~\$142.00~~ Original price was: \$142.00. **\$128.00** Current price is: \$128.00.



**SKU:** RG4070WB

## **PRODUCT DESCRIPTION**

### **PURPOSE AND BENEFITS**

A general purpose, height adjustable, folding walking frame with front wheels and rear skis for easier handling and manoeuvrability over different surfaces. Made of lightweight aluminium, this frame can help improve balance and stability.

### **FEATURES**

- Aluminium anodised frame
- Height adjustable
- Soft foam or rubber handgrips optional attachment
- Swivel Wheels