




### REGENT COMMODE SLING

**\$348.00 - \$448.00** Price range: \$348.00 through \$448.00



**SKU:** 1001141

## PRODUCT DESCRIPTION

### FEATURES AND BENEFITS

- Excellent back support and comfort.
- Suitable for patients who generally have good trunk and hip control.
- Long divided legs.
- Elasticated neoprene waist strap.
- Easy access for clothing removal during toileting.
- Available in polyester solid only.
- Positioning straps.
- Colour coded lifting straps.
- Made to Australian standards.
- N.A.T.A.\* Tested safe working load is 250kg.

For sizing see dimensions (cm) in the diagram below:

### SIZE

- Small A: 800 B: 900 C: 650 D: 200
- Medium A: 850 B: 900 C: 650 D: 200
- Large A: 950 B: 900 C: 650 D: 200
- XL A: 1050 B: 1050 C: 750 D: 215