




REGENT FULL BODY SLING

\$322.00 - \$424.00 Price range: \$322.00 through \$424.00



SKU: SUPSRF

PRODUCT DESCRIPTION

FEATURES AND BENEFITS

- Excellent back and neck support
- For patients who need full body and neck support
- Available in polyester solid and mesh
- Positioning straps
- Colour coded lifting straps
- Made to Australian standards
- N.A.T.A.* Tested safe working load is 250kg

For sizing see dimensions (cm) in the diagram below:

SIZE

- X-Small A: 1050 B: 350 C: 600
- Small A: 1250 B: 400 C: 800
- Medium A: 1300 B: 500 C: 900
- Large A: 1400 B: 550 C: 950