



REGENT GENERAL PURPOSE SLING WITH HEAD SUPPORT

\$305.00 - \$353.00 Price range: \$305.00 through \$353.00



SKU: SUPSGH

PRODUCT DESCRIPTION

FEATURES AND BENEFITS

- For patients who need full body and neck support
- Divided legs can be crossed in the centre, uncrossed or positioned under thighs.
- Available in solid polyester and mesh variants.
- Positioning straps.
- Colour coded lifting straps.
- Made to Australian standards.
- N.A.T.A.* Tested safe working load is 250kg.

For sizing see dimensions (cm) in the diagram below:

SIZE

- X-Small: A: 1100 B: 500 C: 600
- Small A: 1200 B: 550 C: 700
- Medium A: 1400 B: 600 C: 750
- Large A: 1500 B: 650 C: 850