



REGENT GENERAL PURPOSE SLING

\$254.00 - \$322.00 Price range: \$254.00 through \$322.00



SKU: SUPSGP

PRODUCT DESCRIPTION

FEATURES AND BENEFITS

- Divided legs can be crossed in the centre, uncrossed or positioned under thighs
- Suitable for patients with good head control
- Available in polyester solid and mesh variants.
- Positioning straps
- Colour coded lifting straps
- Made to Australian standards
- N.A.T.A.* Tested safe working load is 250kg

For sizing see dimensions (cm) in the diagram below:

SIZE

- X-Small: A: 750 B: 900 C: 500 D: 300
- Small A: 800 B: 1150 C: 550 D: 400
- Medium A: 900 B: 1300 C: 600 D: 450
- Large A: 950 B: 1350 C: 650 D: 500