



### REGENT HALF BODY SLING

**\$263.00 - \$314.00** Price range: \$263.00 through \$314.00



**SKU:** N/A

## PRODUCT DESCRIPTION

### FEATURES AND BENEFITS

- Suitable for patients with good head control
- Quick and easy to fit
- Available in polyester solid and mesh
- Positioning straps
- Colour coded lifting straps
- Made to Australian standards

N.A.T.A.\* Tested safe working load is 250kg

For sizing see dimensions (cm) in the diagram below:

### SIZE

- X-Small A: 550 B: 800 C: 350
- Small A: 750 B: 900 C: 400
- Medium A: 800 B: 1050 C: 500
- Large A: 850 B: 1100 C: 550