

REGENT SLING EXTENDER BELTS

\$64.00 - \$77.00 Price range: \$64.00 through \$77.00



SKU: SUPSEXT-H



PRODUCT DESCRIPTION

Purpose & Benefits of Regent Sling Extender Belts

Regent Sling Extender Belts are designed to provide additional length to patient lifting slings, enhancing their versatility and functionality. These belts are essential in healthcare settings, allowing for safer and more comfortable patient transfers. They help accommodate different patient sizes and lifting needs, ensuring that slings can be adjusted for optimal positioning and support.

Features

- **Adjustable Length:** Available in regular (40 cm) and long (80 cm) sizes to suit various patient and equipment sizes.
- **Durable Neoprene Material:** Ensures strength and longevity, suitable for frequent use in medical environments.
- **Easy Integration:** Compatible with Commode, Deluxe Commode, and Toileting Slings, facilitating a seamless addition to existing sling setups.

Specifications

- **Material:** Neoprene, known for its durability and flexibility.
- **Sizes:** Two lengths available, 40 cm and 80 cm, both 16 cm wide.

These belts are particularly useful in settings that require adaptable solutions for patient handling, providing essential extensions to enhance the functionality of lifting slings.