



REGENT TOILET / STANDING SLING

\$334.00 - \$415.00 Price range: \$334.00 through \$415.00



SKU: SUPSTS

PRODUCT DESCRIPTION

FEATURES AND BENEFITS

- Soft padding throughout sling
- Allows for safe and comfortable chair to toilet transfer
- Suitable for patients who generally have good trunk and hip control
- Allows patient to be easily brought to a standing position
- Elasticated neoprene waist strap
- Easy access for clothing removal during toileting
- Available in solid polyester only
- Positioning straps
- Colour coded lifting straps
- Made to Australian standards
- N.A.T.A.* Tested safe working load is 250kg

SIZE

- Small/Medium A: 965 B: 540 C: 380
- Large/X-Large A: 1080 B: 620 C: 380