




RESMED AIRFIT N20 CUSHION

~~\$88.00~~ Original price was: ~~\$88.00~~. ~~\$54.00~~ Current price is: \$54.00. Sale Available Online Only! While Stocks Last



SKU: 63552-H

PRODUCT DESCRIPTION

Purpose & Benefits

The **ResMed AirFit N20 Cushion** is a replacement nasal cushion designed for use with the AirFit N20 and AirTouch N20 mask systems. Its InfinitySeal™ silicone technology adapts to facial contours and movement to maintain a consistent seal throughout the night, helping reduce leaks and increase comfort. It enables users to benefit from reliable CPAP therapy across different pressures and head movements.

Features

- InfinitySeal™ silicone cushion that adjusts to facial movements for a resilient seal
- Modular cushion design: Small, Medium and Large sizes snap into the same mask frame
- Soft silicone material for comfort and durability with frequent use

Specifications

- **Available Sizes:** Small, Medium, Large
- **Nose Width Guide:** Small ~ 23–37 mm, Large ~ 45–55 mm