



RESMED AIRFIT N20 NASAL MASK FOR HER



SKU: 63520

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

This is a compact nasal mask for CPAP sleep apnea therapy, designed to fit every woman. Comfort and freedom through innovative technology. Easy to use, unobtrusive, suitable for all types of sleep therapy.

FEATURES

- ResMed's 'For Her' range of masks and devices have been specially adapted to suit women's facial features. Experience the same performance as our standard range; enjoy the luxury of a design tailored for women.
- Thanks to its innovative silicone cushion, AirFit N20 for Her is designed to deliver a robust seal for women of all facial shapes, sizes and ethnic backgrounds.
- our advanced InfinitySeal™ silicone cushion adjusts to your movements for comfort and security throughout the night.
- with a mask this simple (magnetic clips, quick release elbow), fitting it and keeping it clean are child's play.
- the under-the-eye frame design provides stability, minimal facial contact and unrestricted vision designed to help you see more, do more and feel great.
- a short flexible tube interconnects the elbow to the machine tube, allowing you to move more easily without encountering tube drag or leakage.
- the flexible frame is lined with plush fabric so it's soft on your skin and remains comfortable throughout the night.
- new technologies, extensive user testing, thoughtful design: we did the hard work so you can enjoy a simple, hassle-free experience.

- It delivers performance with high pressure, Bi-level and NIV therapies as well as standard CPAP.

Warning

Keep the mask magnets at a safe distance of at least 6 inches (150 mm) away from implants or medical devices that may be adversely affected by magnetic interference. This warning applies to you or anyone in close physical contact with your mask. The magnets are in the frame and lower headgear clips, with a magnetic field strength of up to 400mT. When worn, they connect to secure the mask but may inadvertently detach while asleep. Implants/medical devices, including those listed within contraindications, may be adversely affected if they change function under external magnetic fields or contain ferromagnetic materials that attract/repel to magnetic fields (some metallic implants, e.g., contact lenses with metal, dental implants, metallic cranial plates, screws, burr hole covers, and bone substitute devices). Consult your physician and manufacturer of your implant / other medical device for information on the potential adverse effects of magnetic fields.

[Click to view all ResMed affected masks and magnet locations](#)

[Click to open Patient Letter - Urgent Product Defect Correction ResMed Masks with Magnets](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au