


RESMED AIRFIT N20 NASAL MASK



SKU: 6372-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Resmed AirFit N20 is a compact nasal mask for CPAP sleep apnea therapy, designed to fit everyone. Comfort and freedom through innovative technology. Easy to use, unobtrusive, suitable for all types of sleep therapy.

Designed to fit any face shape or size, the AirFit N20 moves comfortably with you through the night, providing a secure fit and comfortable seal tailored to your therapy needs.

KEY FEATURES

- Thanks to its innovative silicone cushion, AirFit N20 is designed to deliver a robust seal for all facial shapes, sizes and ethnic backgrounds.
- our advanced InfinitySeal™ silicone cushion adjusts to your movements for comfort and security throughout the night.
- with a ResMed mask this simple (magnetic clips, quick-release elbow), fitting it and keeping it clean are easy as 1, 2, 3.
- The under-the-eye frame design provides stability, minimal facial contact and unrestricted vision designed to help you see more, do more and feel great.
- a short flexible tube interconnects the elbow to the machine tube, allowing you to move more easily without encountering tube drag or leakage.
- the flexible frame is lined with plush fabric so it's soft on your skin and remains comfortable throughout the night.
- extensive user testing, ground-breaking technologies, thoughtful design: we did the hard work so you can enjoy a hassle-free and simple experience.

- AirFit N20 delivers performance with high pressure, Bi-level and NIV therapies as well as standard CPAP.

AIRFIT N20 GUIDES

- [Mask fitting video guide for the AirFit N20](#)

Buy the Resmed AirFit N20 online at Medimart with AfterPay or NDIS Payment method for your convenience.

Warning

Keep the mask magnets at a safe distance of at least 6 inches (150 mm) away from implants or medical devices that may be adversely affected by magnetic interference. This warning applies to you or anyone in close physical contact with your mask. The magnets are in the frame and lower headgear clips, with a magnetic field strength of up to 400mT. When worn, they connect to secure the mask but may inadvertently detach while asleep. Implants/medical devices, including those listed within contraindications, may be adversely affected if they change function under external magnetic fields or contain ferromagnetic materials that attract/repel to magnetic fields (some metallic implants, e.g., contact lenses with metal, dental implants, metallic cranial plates, screws, burr hole covers, and bone substitute devices). Consult your physician and manufacturer of your implant / other medical device for information on the potential adverse effects of magnetic fields.

[Click to view all ResMed affected masks and magnet locations](#)

[Click to open Patient Letter - Urgent Product Defect Correction ResMed Masks with Magnets](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au