




## RESMED AIRFIT P10 NASAL PILLOW SYSTEM



**SKU:** 62904

## PRODUCT DESCRIPTION

### RESMED AIRFIT P10 PURPOSE & BENEFITS

Resmed AirFit P10 nasal pillows mask is our quietest mask yet! Bringing you the latest technology in sleep apnea treatment, the AirFit P10 retains the proven performance and reliability of the popular Swift™ FX mask series but is approximately 50% quieter and 50% lighter.

Providing a stable seal and a high level of comfort, the AirFit P10 is also easy to use, with just three pieces for quick assembly and cleaning. The mask's innovative QuickFit™ headgear is unobtrusive yet snug, giving you a real sense of facial freedom while you're wearing it. That means you can easily shift to your favourite position without interrupting your nightly sleep apnea therapy.

### KEY FEATURES

- Easy to use. The AirFit P10's minimalist three-piece design makes assembly, fitting and cleaning a breeze.
- Easy to assemble. Key-clip nasal pillows feature colour-coded sizing and left/right-orientation for fast and easy fitting. The mask easily slips on and off without any extra steps.
- Comfortable. The mask is designed to seal on contact and lets you move around comfortably without compromising the seal
- Suitable for higher pressures

Nasal pillows give minimum facial coverage and are lightweight. Patients additionally enjoy the fact these masks don't have a lot of matter intruding on their faces. A good quality Nasal pillow reduces air leakage because they deliver air into the patient's nasal passage.

CPAP masks don't get much smaller and simpler than AirFit P10 – weighing only 1.6 oz. and made up of only three parts! With a nasal pillows cushion and slip-on headgear, it's easy to use and covers very little of your face. The new QuickFit headgear clip makes it easy to adjust at home, and for a more comfortable night's sleep, it features QuietAir woven-mesh vents that gently disperse exhaled air in multiple directions away from you and your partner.

For more products like this AirFit P10 nasal pillow Check out our [CPAP Machines](#).

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)