


RESMED AIRSENSE 10 CLIMATELINEAIR HEATED TUBE



SKU: 37296



PRODUCT DESCRIPTION

The heated tube for the [Airsense](#)™ 10 and AirCurve™ 10 line by ResMed, the ClimateLineAir™ Heated Tube, provides comfort and optimal humidity by monitoring the tube temperature at intervals between the mask and the machine which allow the humidifier to compensate as the ambient conditions change. The innovative swivel connection at the back of the machine allows for comfortable movement through the night.

The ClimateLineAir™ Heated Tube for AirSense™ and AirCurve™ heated hoses keep warm, humidified air all the way to the mask while at the same time avoiding rainout.

Once you connect the ClimateLineAir™ Heated Tube to an AirSense or AirCurve 10 machine, the Auto Climate Control will engage automatically. The machine will monitor ambient humidity levels and then alter the moisture and tube temperature settings to ensure your night is as comfortable as possible.

The 360 degree swivel feature of the ClimateLineAir™ allows freedom of movement while keeping the electrical connection secure.

- This tubing is ONLY compatible with the ResMed Airsense™ and AirCurve™ 10.

For further information go to **[medimart.com.au](https://www.medimart.com.au)** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au