



RESMED CLIMATELINE HEATED TUBE FOR S9



SKU: 36995

PRODUCT DESCRIPTION

The ClimateLine tube enhances the H5i Humidifier by controlling the air temperature as it passes through the tube. This heated tube is compatible with the ResMed S9 range.

By keeping the temperature in the tube steady, the Climate Control system automatically adjusts to changes in ambient conditions to maintain a comfortable humidity level. Rainout is avoided without needing to reduce the humidity delivered. A temperature sensor built into the mask end of the tube is automatically controlled by the Climate Control algorithm to deliver a consistent, comfortable temperature at the mask.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au