



### RESMED HUMIDX F20 3-PACK

~~\$43.00~~ Original price was: ~~\$43.00~~. **\$37.00** Current price is: **\$37.00**. Sale Available Online Only! While Stocks Last



**SKU:** 38012

## PRODUCT DESCRIPTION

### Purpose & Benefits

The **ResMed HumidX F20 3-Pack** provides waterless humidification for AirMini users with F20 full face masks. Designed to retain heat and moisture, it delivers effective humidification without the need for a traditional water chamber, making it ideal for travel and convenience.

### Features

- Waterless humidification system designed for ResMed AirMini with F20 masks
- Captures heat and moisture to provide comfortable CPAP therapy
- Compact design for travel and daily use

### Specifications

- **Pack Size:** 3 HumidX F20 cartridges
- **Compatibility:** ResMed AirMini with AirFit F20 and AirTouch F20 masks