



REVOLUTION SWIVEL CHAIR



SKU: REVO-H

PRODUCT DESCRIPTION

REVOLUTION SWIVEL CHAIR PURPOSE & BENEFITS

The Revolution Swivel Chair is a unique and innovative chair that allows easy access to the seat and to a table edge without moving its legs. The chair incorporates a patented "swivel-slide" mechanism that allows the seat to turn and slide, while the legs remain stationary. This makes it easy for people of all abilities to sit down and get up from the chair, even if they are wheelchair-bound or have limited mobility.

The Revolution Swivel Chair's seat and back are made of commercial-grade vinyl, which is easy to clean and disinfect. The chair is also sturdy and durable, making it a great choice for a variety of settings, including dining rooms, rehabilitation facilities, and in-home care.

FEATURES

- The Revolution Chair incorporates a unique, patented "swivel-slide" technology
- Allows easy access to the chair and to a table edge without moving its legs
- Standard is height adjustable, the bariatric is fixed.
- While the chair legs are stationary, the seat turns to receive the user and returns to face the table
- It then glides forward to position the user at a comfortable distance from the table edge
- The lock/release mechanism allows the seat to gently slide back from the table so it can swivel left or right for easy exit

If you are looking for a comfortable and easy-to-use chair that is perfect for people of all abilities, then the Revolution Swivel Chair is the perfect choice for you.

SPECIFICATIONS

Standard

- Max Seat Height: 530mm
- Min Seat Height: 450mm
- Seat Width: 45cm
- Seat Depth: 48cm
- Seat to Floor Height: 480mm
- SWL: 125kg

Bariatric

- Seat Width: 660mm
- Seat Depth: 450mm
- Overall Width: 790mm
- Overall Depth: 570mm
- Overall Height: 850mm
- Seat to Floor Height: 480mm
- SWL: 250kg

Height Adjustable

- Total Width: 505mm
- Total Depth: 560mm
- Seat Depth: 485mm
- Seat Width: 450mm
- Seat to Floor Height: 450mm - 530mm
- SWL: 120kg

OPTIONAL EXTRAS

- 50mm Wider
- [Back Extension Brackets](#)
- [Easy Lever](#)
- Front Locking Setup
- [Non-slip Rubber feet](#)
- High back

Please [contact us](#) for any customisations or optional extras.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au