



REVOLUTION TIMBER SWIVEL CHAIR

\$1,795.00



SKU: REVO-T-H

PRODUCT DESCRIPTION

REVOLUTION SWIVEL CHAIR PURPOSE & BENEFITS

The Revolution Swivel Chair is a unique and innovative chair that allows easy access to the seat and to a table edge without moving its legs. The chair incorporates a patented "swivel-slide" mechanism that allows the seat to turn and slide, while the legs remain stationary. This makes it easy for people of all abilities to sit down and get up from the chair, even if they are wheelchair-bound or have limited mobility.

The Revolution Swivel Chair's seat and back are made of commercial-grade vinyl, which is easy to clean and disinfect. The chair is also sturdy and durable, making it a great choice for a variety of settings, including dining rooms, rehabilitation facilities, and in-home care.

FEATURES

- The Revolution Chair incorporates a unique, patented "swivel-slide" technology
- Allows easy access to the chair and to a table edge without moving its legs
- Standard is height adjustable.
- While the chair legs are stationary, the seat turns to receive the user and returns to face the table
- It then glides forward to position the user at a comfortable distance from the table edge
- The lock/release mechanism allows the seat to gently slide back from the table so it can swivel left or right for easy exit

If you are looking for a comfortable and easy-to-use chair that is perfect for people of all abilities, then the Revolution Swivel Chair is the perfect choice for you.

SPECIFICATIONS

Standard

- Max Seat Height: 530 mm
- Min Seat Height: 450 mm
- Seat Width: 45 cm
- Seat Depth: 48 cm
- SWL: 125 kg

OPTIONAL EXTRAS

- 50mm Wider
- Back Extension Brackets
- Easy Lever
- Front Locking Setup
- Non-slip Rubber feet
- High back

Please [contact us](#) for any customisations or optional extras.