

ROHO ADAPTOR PAD

\$375.00



SKU: ADAPTOR



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The self-contained air cells in the ADAPTOR Pad can be used as an interface on virtually any flat or contoured surface to protect skin and soft tissues. Especially useful when cut to fit various surface shapes such as sides of wheelchairs, trays, head and footrests, inside of helmets, or for any surfaces that may need protection.

FEATURES

- Absorbs shock
- Protects against shear and friction
- Easily cleaned and disinfected
- Durable
- Closed air cells, so product maintains efficacy even if a single cell is cut or punctured
- Made of flame-resistant, black neoprene (not made with natural rubber latex)

Protection

ROHO is the worldwide leader of seating solutions to prevent and treat pressure injuries, including deep tissue injury, while providing skin and soft tissue protection.

Clinical Evidence

ROHO has been in over 90 scientific and clinical studies to reinforce and verify that a ROHO DRY FLOATATION cushion heals, treats and prevents pressure injuries, including deep tissue injury.

Adjustability

ROHO DRY FLOATATION cushions, unlike foam products, provide the ability to constantly adjust and conform with an individual's body, now, and as he or she changes throughout the months and years.

Adaptability

ROHO DRY FLOATATION cushions support the shifts and movements that individuals make throughout each day.

SPECIFICATIONS

- Width: 9.75" (250mm)
- Depth: 13.5" (340mm)
- Height: 1" (20mm)
- Product Weight: 0.2kg
- User Weight Limit: No limit

[View the ROHO Cushions Selection Guide](#)