


## ROHO ADAPTOR PAD



**SKU:** ADAPTOR

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The self-contained air cells in the ADAPTOR Pad can be used as an interface on virtually any flat or contoured surface to protect skin and soft tissues. Especially useful when cut to fit various surface shapes such as sides of wheelchairs, trays, head and footrests, inside of helmets, or for any surfaces that may need protection.

### FEATURES

- Absorbs shock
- Protects against shear and friction
- Easily cleaned and disinfected
- Durable
- Closed air cells, so product maintains efficacy even if a single cell is cut or punctured
- Made of flame-resistant, black neoprene (not made with natural rubber latex)

### Protection

ROHO is the worldwide leader of seating solutions to prevent and treat pressure injuries, including deep tissue injury, while providing skin and soft tissue protection.

### Clinical Evidence

ROHO has been in over 90 scientific and clinical studies to reinforce and verify that a ROHO DRY FLOATATION cushion heals, treats and prevents pressure injuries, including deep tissue injury.

### Adjustability

ROHO DRY FLOATATION cushions, unlike foam products, provide the ability to constantly adjust and conform with an individual's body, now, and as he or she changes throughout the months and years.

**Adaptability**

ROHO DRY FLOATATION cushions support the shifts and movements that individuals make throughout each day.

**SPECIFICATIONS**

- Width: 9.75" (250mm)
- Depth: 13.5" (340mm)
- Height: 1" (20mm)
- Product Weight: 0.2kg
- User Weight Limit: No limit

[View the ROHO Cushions Selection Guide](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

**SA****Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

**Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

**Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

**VIC****Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

**Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

**NT**

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)