


ROHO AIRLITE CUSHION

\$390.00 - \$410.00 Price range: \$390.00 through \$410.00



SKU: AL-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The AirLITE cushion is a non-adjustable, maintenance-free foam cushion with a ROHO air insert that conforms to your body while also preventing you from bottoming out in your cushion. The gentle contour of the leg trough helps you maintain the right position and posture, and lets you feel more secure.

FEATURES

- Designed for clients with tactile sensation who are at low risk of developing skin breakdown
- Contoured foam base with an internal sealed air cell provides comfort and support for basic positioning and pressure management needs.
- **“Bottom-Out” Protection:** The non-adjustable ROHO AIR FLOATATION™ air insert built into the cushion under the ischial has been enhanced to provide better bottom-out protection.
- **Softer Foam:** Softer foam allows the cushion to conform to the individual’s physiology and provide a more comfortable seating surface.
- **Durability and Simplicity:** The 24 month warranty shows our confidence in the quality and durability of the product. The AirLITE is an off the shelf, ready to use product. No adjustment is required.
- **Heavy Duty Cover:** The cover design, with full zipper coverage, provides enhanced fluid resistance to protect the foam. At the same time, the cover is 100% PVC free, reinforcing ROHO’s commitment to more environmentally friendly products.
- All models include: 2 way stretch cover, Hand pump & Repair kit, Operating instructions.

SPECIFICATIONS

Product	Total Height	Safe Working Load	Pressure Care Rating	Material
ROHO AirLite 2 - 3.75 in		113 kg	Low Risk	Air Cells

[View the ROHO Cushions Selection Guide](#)