

ROHO COMMUNE CUSHION LOW PROFILE

\$1,195.00 - \$1,310.00 Price range: \$1,195.00 through \$1,310.00



SKU: SEATCMD88LP-H

PRODUCT DESCRIPTION

Purpose and benefits

The ROHO commode cushion low profile is designed to deliver effective pressure relief and skin protection for individuals using shower chairs or commodes. Built with ROHO's DRY FLOATATION® technology, its air-filled cells contour to the user's shape, helping to evenly distribute pressure and reduce the risk of skin breakdown during hygiene routines.

The low-profile design supports safe transfers and promotes comfort while seated, making it ideal for those at risk of pressure injuries or requiring prolonged use of commode seating.

Features

- 5.7 cm low-profile air cells enhance stability and simplify transfers
- ROHO DRY FLOATATION® technology offers uniform pressure distribution
- Made from flame-resistant, black neoprene (not made with natural rubber latex)
- Waterproof, anti-static, and easy to clean
- Lightweight and portable for home, facility, or travel use
- Compatible with standard commodes and shower chairs

- Supplied with hand pump, repair kit, and user instructions

Specifications

- Cell height: 5.7 cm
- Available sizes (width x depth): 37.5 cm x 37.5 cm, 42 cm x 42 cm, 46.5 cm x 46.5 cm
- Material: Flame-resistant black neoprene
- Weight: Approx. 0.7 kg
- Latex: Not made with natural rubber latex
- Weight capacity: No limit when properly sized
- Warranty: 24-month limited

[View the ROHO Cushions Selection Guide](#)