



ROHO CONTOUR BASE

\$160.00



SKU: BASE1313-H

PRODUCT DESCRIPTION

Purpose and benefits

The ROHO contour base is a lightweight, moulded insert designed to provide a stable seating foundation for individuals using sling-seat wheelchairs. It helps convert the hammock-style seat into a more supportive and symmetrical surface, improving posture and comfort when used with a ROHO cushion.

By promoting forward pelvic stability and preventing femoral rotation, the contour base enhances seated alignment and reduces the risk of postural issues. It sits discreetly beneath the cushion within the standard ROHO cushion cover, ensuring minimal impact on overall cushion performance while adding structural support.

Features

- Provides a rigid base for sling-seat wheelchairs
- Promotes symmetrical seating and improved pelvic alignment
- Helps prevent femoral rotation and anterior pelvic tilt
- Sits inside the ROHO cushion cover for seamless integration
- Made from durable, lightweight moulded plastic

- Easy to insert and remove as needed
- Enhances stability without altering cushion pressure-relief function

Specifications

- Material: Moulded plastic
- Compatibility: Fits within ROHO cushion covers for sling-seat wheelchairs
- Weight: Approx. 0.7 kg
- Warranty: 12-month limited

[View the ROHO Cushions Selection Guide](#)