



ROHO CONTOUR BASE



SKU: BASE1313-H

PRODUCT DESCRIPTION

Purpose and benefits

The ROHO contour base is a lightweight, moulded insert designed to provide a stable seating foundation for individuals using sling-seat wheelchairs. It helps convert the hammock-style seat into a more supportive and symmetrical surface, improving posture and comfort when used with a ROHO cushion.

By promoting forward pelvic stability and preventing femoral rotation, the contour base enhances seated alignment and reduces the risk of postural issues. It sits discreetly beneath the cushion within the standard ROHO cushion cover, ensuring minimal impact on overall cushion performance while adding structural support.

Features

- Provides a rigid base for sling-seat wheelchairs
- Promotes symmetrical seating and improved pelvic alignment
- Helps prevent femoral rotation and anterior pelvic tilt
- Sits inside the ROHO cushion cover for seamless integration
- Made from durable, lightweight moulded plastic

- Easy to insert and remove as needed
- Enhances stability without altering cushion pressure-relief function

Specifications

- Material: Moulded plastic
- Compatibility: Fits within ROHO cushion covers for sling-seat wheelchairs
- Weight: Approx. 0.7 kg
- Warranty: 12-month limited

[View the ROHO Cushions Selection Guide](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au