

ROHO CONTOUR SELECT CUSHION



SKU: CS-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the revolutionary Roho Contour Select Cushion, meticulously crafted to enhance your comfort, posture, and skin protection. Engineered with ingenious precision, this cushion boasts multiple cell heights that cater to your unique positioning needs. Whether you're seeking optimal posture support or superior skin protection, the Contour Select has got you covered.

One of its standout features is the exclusive ISOFLO Memory Control, empowering you to customize your seating experience like never before. By simply locking air into the four distinct compartments, this innovative technology delivers exceptional stability and precise positioning, effectively minimizing any unwanted side-to-side or front-to-back motion. Say goodbye to discomfort and hello to a newfound sense of stability and support with the Roho Contour Select Cushion.

FEATURES

- A unique combination of high (4"), mid (3") and low(2") profile cells in the one cushion centre the user and position the pelvis back in the wheelchair providing greater lateral stability through the channelling of the legs
- ISOFLO Memory Control locks air into compartments to help maintain better postural integrity in addition to pressure reduction
- Interconnected cells help disburse pressure and reduce shear risk promoting pressure management
- Appropriate for clients with a high risk of ischemic ulcers and need for greater positioning support
- Soft neoprene construction is flame-retardant and durable

SPECIFICATIONS

Height	Approx. 2.25" - 4.25" 5.5cm - 10.5 cm (range due to varying heights of cells)
Weight limit	No weight limit

[View the ROHO Cushions Selection Guide](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au