

ROHO HYBRID ELITE 1R CUSHION SENSOR READY



SKU: 1RHE1414-SR-H

PRODUCT DESCRIPTION

Purpose and benefits

The ROHO hybrid elite 1r cushion is engineered to deliver both advanced skin protection and optimal postural support for wheelchair users. This cushion combines ROHO's DRY FLOATATION® air cell technology with a contoured foam base, making it ideal for users who require enhanced stability without compromising pressure relief.

Its single-compartment air insert molds to the user's body, evenly redistributing pressure and reducing shear and friction forces. The contoured foam base promotes pelvic alignment and provides targeted support for the thighs, making it suitable for both symmetrical and asymmetrical postures.

Features

- DRY FLOATATION® air cells adjust to the user for optimal pressure redistribution
- Contoured foam base supports pelvic alignment and thigh positioning
- Single-valve air compartment for simple and precise inflation
- Reduces shear and friction for improved skin integrity
- Compatible with optional Smart Check device for inflation monitoring

- Includes hand pump, repair kit, and cover
- Washable cover and wipe-clean foam base

Specifications

- Available sizes (width × depth): 34.5 cm × 37.5 cm up to 60 cm × 52.5 cm
- Cushion height: approx. 11 cm (unloaded)
- Cushion weight: approx. 1.4 kg (varies by size)
- Maximum user weight: 225 kg
- Materials: Flame-retardant neoprene air insert, high-density foam base
- Warranty: 24 months (cushion), 6 months (cover)

[View the ROHO Cushions Selection Guide](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au