

ROHO HYBRID SELECT CUSHION



SKU: HS-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Hybrid Select is ROHO's most advanced hybrid cushion that provides you with superior positioning and stability, plus the ability to offload by simply removing the IT air insert. You can easily find the right position for optimal immersion and skin protection by utilising the exclusive ISO FLO Memory Control. This technology allows you to set the amount of air in the three separate ROHO compartments—the removable IT insert and the dual-zone overlay.

The Hybrid Select also uses the positioning capabilities of a contoured, high resiliency foam base with deep leg troughs to offer you support, lateral stability, and provide relief to the coccyx.

FEATURES

- Stabilises and holds the position of your pelvis and thighs
- Unique ROHO air cell design specifically crafted to work with the deeper curves of the contoured foam base
- Delivers skin and soft tissue protection, with optional offloading under the ischial tuberosities (ITs)
- Easy to set up

SPECIFICATIONS

Cover type	Stretch-Air
-------------------	-------------

Product weight	3.2lbs (1.45 kg)*
Weight limit	500 lbs (227 kg)

[View the ROHO Cushions Selection Guide](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au