


## ROHO HYBRID SELECT CUSHION

**\$1,040.00 - \$1,125.00** Price range: \$1,040.00 through \$1,125.00



**SKU:** HS-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Hybrid Select is ROHO's most advanced hybrid cushion that provides you with superior positioning and stability, plus the ability to offload by simply removing the IT air insert. You can easily find the right position for optimal immersion and skin protection by utilising the exclusive ISO FLO Memory Control. This technology allows you to set the amount of air in the three separate ROHO compartments—the removable IT insert and the dual-zone overlay.

The Hybrid Select also uses the positioning capabilities of a contoured, high resiliency foam base with deep leg troughs to offer you support, lateral stability, and provide relief to the coccyx.

### FEATURES

- Stabilises and holds the position of your pelvis and thighs
- Unique ROHO air cell design specifically crafted to work with the deeper curves of the contoured foam base
- Delivers skin and soft tissue protection, with optional offloading under the ischial tuberosities (ITs)
- Easy to set up

### SPECIFICATIONS

<b>Cover type</b>	Stretch-Air
-------------------	-------------

<b>Product weight</b>	3.2lbs (1.45 kg)*
<b>Weight limit</b>	500 lbs (227 kg)

[View the ROHO Cushions Selection Guide](#)