




ROHO LTV LOVE THAT VALVE SEAT CUSHION

\$320.00 - \$410.00 Price range: \$320.00 through \$410.00



SKU: LTV-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Roho LTV Seat Cushion allows you to customize the cushion to suit your specific needs. It is a non-medical cushion for anyone who sits for extended periods of time.

The LTV Cushion is constructed of individual, interconnected air cells that allow air to slowly transfer from each chamber and distributes pressure evenly through the entire cushion. By distributing the weight, the pressure is reduced and comfort is maximised.

FEATURES

- For everyday comfort and pressure care
- Adjustable airflow
- Choice of 2 covers, black fabric or black Leatherette

SPECIFICATIONS

- Width: 18" / 45.5 cm
- Depth: 16" / 40.5 cm
- Height: 1" to 2" / 2.5 cm to 5.5 cm

[View the ROHO Cushions Selection Guide](#)