


ROHO MINI-MAX CUSHION



SKU: MINI-H



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the ROHO MINI-MAX Cushion - designed for the active individual seeking superior comfort and support. This remarkable wheelchair cushion offers outstanding shock absorption, skin protection, and stability. Engineered with advanced technology, it ensures an unmatched experience for those with active lifestyles.

FEATURES

- Minimizes shear and friction.
- Provides support and comfort in the scapular sacral, or lumbar region.
- Versatile use. Excellent for sports, bathing, and travel.
- Made of non-permeable, non-static and flame-resistant, black neoprene (not made with natural rubber latex).
- Carry it anywhere.

SPECIFICATIONS

Material	Air Cells
Safe Working Load	Unlimited
Total Height	32 mm

[View the ROHO Cushions Selection Guide](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au