

ROHO PACK-IT CUSHION



SKU: PACKIT

PRODUCT DESCRIPTION

Purpose and Benefits

The ROHO Pack-It Cushion is designed to provide targeted, short-term pressure relief for active users. Compact and lightweight, it is ideal for use during travel, sports, bathing, or as support in scapular, sacral, or lumbar regions.

Its air-cell design promotes pressure redistribution and shock absorption, helping prevent skin breakdown. Whether used on a wheelchair, stadium seat, or standard chair, this cushion delivers versatile comfort in a range of everyday settings.

Features

- Lightweight and compact for easy portability
- Air-cell technology provides pressure redistribution and shock absorption
- Ideal for use as scapular, sacral, or lumbar support
- Made from flame-resistant, black neoprene (not made with natural rubber latex)
- Includes two-way stretch cover, hand inflation pump, repair kit, and instructions
- Easy to clean with soap and water

- Suitable for active use during travel or athletic activity

Specifications

- Dimensions: 41.5 cm x 24 cm x 3 cm
- Weight: Approx. 700 g
- Material: Flame-resistant neoprene
- Latex: Not made with natural rubber latex

[View the ROHO Cushions Selection Guide](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au