

ROHO PACK-IT CUSHION



\$710.00 Original price was: \$710.00.\$639.00 Current price is: \$639.00.



SKU: PACKIT

PRODUCT DESCRIPTION

Purpose and Benefits

The ROHO Pack-It Cushion is designed to provide targeted, short-term pressure relief for active users. Compact and lightweight, it is ideal for use during travel, sports, bathing, or as support in scapular, sacral, or lumbar regions.

Its air-cell design promotes pressure redistribution and shock absorption, helping prevent skin breakdown. Whether used on a wheelchair, stadium seat, or standard chair, this cushion delivers versatile comfort in a range of everyday settings.

Features

- Lightweight and compact for easy portability
- Air-cell technology provides pressure redistribution and shock absorption
- Ideal for use as scapular, sacral, or lumbar support
- Made from flame-resistant, black neoprene (not made with natural rubber latex)
- Includes two-way stretch cover, hand inflation pump, repair kit, and instructions
- Easy to clean with soap and water



• Suitable for active use during travel or athletic activity

Specifications

• Dimensions: 41.5 cm x 24 cm x 3 cm

• Weight: Approx. 700 g

• Material: Flame-resistant neoprene

• Latex: Not made with natural rubber latex

View the ROHO Cushions Selection Guide