

ROHO PACK-IT CUSHION

\$710.00



SKU: PACKIT



PRODUCT DESCRIPTION

Purpose and Benefits

The ROHO Pack-It Cushion is designed to provide targeted, short-term pressure relief for active users. Compact and lightweight, it is ideal for use during travel, sports, bathing, or as support in scapular, sacral, or lumbar regions.

Its air-cell design promotes pressure redistribution and shock absorption, helping prevent skin breakdown. Whether used on a wheelchair, stadium seat, or standard chair, this cushion delivers versatile comfort in a range of everyday settings.

Features

- Lightweight and compact for easy portability
- Air-cell technology provides pressure redistribution and shock absorption
- Ideal for use as scapular, sacral, or lumbar support
- Made from flame-resistant, black neoprene (not made with natural rubber latex)
- Includes two-way stretch cover, hand inflation pump, repair kit, and instructions
- Easy to clean with soap and water

- Suitable for active use during travel or athletic activity

Specifications

- Dimensions: 41.5 cm x 24 cm x 3 cm
- Weight: Approx. 700 g
- Material: Flame-resistant neoprene
- Latex: Not made with natural rubber latex

[View the ROHO Cushions Selection Guide](#)