

ROHO SINGLE COMPARTMENT CUSHION



SKU: 1R88LPC-H



PRODUCT DESCRIPTION

Purpose and benefits

The ROHO single compartment cushion is a premium pressure relief solution designed to enhance comfort, skin protection, and posture support for wheelchair users. Its air-filled, interconnected cell system automatically adapts to the user's body contours, distributing pressure evenly and promoting healthy blood circulation.

Ideal for individuals at risk of pressure injuries or those seated for extended periods, this cushion offers both therapeutic and preventative care. It helps reduce shear and friction, supports stability, and provides the flexibility needed for daily use in various settings.

Features

- Adjustable air-filled cells provide personalised pressure redistribution
- Made from flame-retardant neoprene rubber for durability and hygiene
- Lightweight and portable, ideal for active wheelchair users
- Reduces risk of pressure sores by promoting airflow and reducing heat build-up
- Easy to inflate, adjust, and maintain with included hand pump

- Compatible with a wide variety of wheelchair models
- Supplied with repair kit and instruction manual

Specifications

- **Width:** 11.5 - 23.5"
- **Depth:** 11.5 - 20"
- **Height:** 105mm without load
- **Weight:** Approximately 1.5kg
- **Warranty:** 24 months (limited)

[View the ROHO Cushions Selection Guide](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au