



ROHO SMART CHECK DEVICE

\$335.00



SKU: R-SMARTCHECK-01

PRODUCT DESCRIPTION

Purpose & Benefits

The ROHO® Smart Check is a personal electronic device designed to help users, carers, and clinicians ensure that their ROHO Smart Series cushions are always set at the optimal inflation level. With just the press of a button, Smart Check provides real-time feedback on whether your cushion is properly inflated—removing the guesswork from pressure care.

Ideal for individuals at risk of pressure injuries, this tool supports long-term skin integrity by helping maintain consistent, personalised cushion settings. It's especially valuable for users who depend on high-precision positioning and pressure redistribution in their seating system.

Features

- Works with ROHO Smart Series cushions (e.g., High/Low Profile Smart Cells)
- Provides immediate feedback on cushion inflation levels
- Supports skin protection by maintaining optimal air volume
- Simple one-touch operation with easy-to-read LED indicators
- Assists carers and clinicians in achieving and maintaining correct inflation

- Stores a user's optimal inflation setting for ongoing monitoring
- Promotes user independence and seating confidence
- Lightweight and portable

Specifications

- Compatibility: ROHO Single Compartment Smart Series Cushions
- Function: Measures and monitors air pressure for optimal inflation
- Indicators: Green (correct), Red (too firm), Yellow (too soft)
- Power: Battery-operated (long-lasting, replaceable)
- Memory Function: Stores user-specific optimal setting

[View the ROHO Cushions Selection Guide](#)