



ROHO TOILET SEAT CUSHION

\$740.00



SKU: SEATTOILET

PRODUCT DESCRIPTION

Purpose and Benefits

The ROHO Toilet Seat Cushion is designed to provide enhanced stability, comfort, and skin protection during toileting. Ideal for use on standard toilet seats and portable bedside commodes, this cushion helps reduce the risk of pressure injury by distributing weight evenly and minimizing tissue breakdown.

Its unique multi-height air cell design positions the user securely while supporting sensitive areas, making it suitable for individuals with high skin integrity needs. The cushion attaches easily and securely using adjustable straps, providing a dependable solution for pressure relief during toileting routines.

Features

- Designed for standard toilets and portable bedside commodes
- Multi-height air cells: 5.5 cm outer cells, 2.5 cm inner cells
- Enhances stability and promotes skin protection
- Easy attachment with secure fastening straps
- Made from flame-resistant, anti-static neoprene (not made with natural rubber latex)

- Easy to clean with soap and water
- Lightweight, portable, and ideal for clinical or home use
- Includes two cushion sections, user instructions, and repair kit

Specifications

- Compatibility: Standard-size toilets and portable commodes
- Cell Height: Outer cells 5.5 cm, inner cells 2.5 cm
- Material: Flame-resistant black neoprene, anti-static
- Latex: Not made with natural rubber latex
- Cleaning: Washable with soap and water
- Inclusions: 2 cushion sections, repair kit, user manual
- Colour: Black

[View the ROHO Cushions Selection Guide](#)