




ROHO TOILET SEAT CUSHION



SKU: SEATTOILET

PRODUCT DESCRIPTION

Purpose and Benefits

The ROHO Toilet Seat Cushion is designed to provide enhanced stability, comfort, and skin protection during toileting. Ideal for use on standard toilet seats and portable bedside commodes, this cushion helps reduce the risk of pressure injury by distributing weight evenly and minimizing tissue breakdown.

Its unique multi-height air cell design positions the user securely while supporting sensitive areas, making it suitable for individuals with high skin integrity needs. The cushion attaches easily and securely using adjustable straps, providing a dependable solution for pressure relief during toileting routines.

Features

- Designed for standard toilets and portable bedside commodes
- Multi-height air cells: 5.5 cm outer cells, 2.5 cm inner cells
- Enhances stability and promotes skin protection
- Easy attachment with secure fastening straps
- Made from flame-resistant, anti-static neoprene (not made with natural rubber latex)

- Easy to clean with soap and water
- Lightweight, portable, and ideal for clinical or home use
- Includes two cushion sections, user instructions, and repair kit

Specifications

- Compatibility: Standard-size toilets and portable commodes
- Cell Height: Outer cells 5.5 cm, inner cells 2.5 cm
- Material: Flame-resistant black neoprene, anti-static
- Latex: Not made with natural rubber latex
- Cleaning: Washable with soap and water
- Inclusions: 2 cushion sections, repair kit, user manual
- Colour: Black

[View the ROHO Cushions Selection Guide](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au