



## ROLLZ FLEX - THE INCOGNITO WALKING FRAME

**\$995.00**



**SKU: 1510-H**

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Rollz Flex is a world first for its unique design and function and is perfect for those who need some support of a walking frame but want an alternative look to a traditional walker. Its large bag can safely carry up to 20kg so is a great companion for any mission. Light, strong, easy to steer and stylish.

### FEATURES

- Walking with large bag - the Flex boasts a large 25 litre shopping bag. The bag's centrally located position on the frame means that even when loaded with up to 20kg of shopping the frame will not tip. Compare this to other frames with the bag or basket located on the front of the frame which can be unsafe. This bag also includes an inner bag which can be easily lifted in and out.
- The highest quality fittings including sustainable, hard-wearing waterproof material on the charcoal bag.
- Comfortable padded and stylish seat using faux leather in black.
- Soft tyres made from a PU/foam mix for natural suspension
- 8inch wheels and short wheel base which make it very manoeuvrable and kind to your back

### ATTACHMENTS

[Click here](#) to view our full range of Walker/Frame attachments.

<https://youtu.be/iyW2G1ooZQ0>

## **SPECIFICATIONS**

- Rollz Flex: 7.7kg
- Rollz Flex bag: 25 litre
- Maximum carrying load: 125 kg
- Minimum and maximum height of user: 1.55- 1.95cm
- Height: 80-100cm
- Length: 60 cm
- Width: 62.3 cm
- Folded width: 33cm
- Seat height: 62 cm
- Seat width: 44 cm
- Diameter wheel: 8 inches